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Vacuum Transfer System for powders and grains

AGIERRE[®]
Engineering & Innovative Solutions

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The Agierre's Pneumatic Transfer Systems are the result of the most advanced technology in the handling of powders.

They are mainly applied both in the productive department (on mixers, dryers, fluid bed coolers, blending systems) and in the packaging units (on packaging machines, tableting press machines, etc).

The functioning is based on the suction (by vacuum) of powders and grains that must be transferred from a container to another one.

The Vacuum Pump (modelled according to the Venturi Principle) creates vacuum and depression inside the Vacuum Conveyor and consequently the product suction. When the product is inside the Vacuum Conveyor, it will be separated from the air by a filter. After that, the product falls down thanks to gravity and centrifugal force. The main functioning is based on the compressed air, so the contamination of products is minimum. The Vacuum allows the handling of the powder with a low velocity which guaranties a constant and consistent handling of powders without shock and de-mixing.

On the Vacuum Transfer System can be installed other pumps such us liquid ring pump or similar.

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The Vacuum Transfer Systems have an innovative design and they can be customized according to the technical requirements of each Customer. According to the technical and processing specification the Vacuum Transfer Systems can have different configurations:

- Basic;
- Overpressure;
- Inertisation;

The connections could be with clamp, wing nuts or flanged (also in the UNI PN16 version).

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AUTOMATIC DISPENSING PLANT
FOR RAW MATERIALS

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Control Panel:
- pneumatic
- PLC system

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The Vacuum Conveyors (patent pending) are manufactured in compliance with the "cGMP guidelines", the UNI EN ISO norms and the ATEX and Machine Directives. According to the AGIERRE's manufacturing philosophy, all Conveyors can be customized in according to specific needs.

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The Conveyor could be also used in the applications with explosive powders which require, during the discharge phase, an oxygen residual value inferior to 6%. It's possible with the inertisation process.

Sistemi di trasporto in vuoto con sistema C.I.P. e S.I.P. Vacuum Transfer System with C.I.P. and S.I.P. systems

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These cleaning systems find their main application in the pharmaceutical and chemical industries, particularly in the sterile production areas where it's needed to preserve the chemical and physical properties of the product and avoid the "cross contamination".

The Vacuum Transfer Systems can be equipped with some accessories that allow to wash, dry and sterilize it through the C.I.P. (cleaning in place) and S.I.P. (steam in place) systems. The C.I.P. system cleans the filter and all internal surfaces of Vacuum Conveyor by rotating nozzles (patent pending system).

The C.I.P. system is based on a mechanical cleaning system that ensure the complete removal of all product residues inside the Vacuum Conveyor. After washing phase all surfaces can be dried with air and/or nitrogen. This Conveyor has an innovative and specific design without any point of product or liquid accumulation.

The Vacuum Conveyors can be managed by a pneumatic control console or by a PLC system.

They are manufactured with stainless steel AISI 316L and are in compliance with the "cGMP guidelines", the UNI EN ISO norms, the Atex and Machine Directives.

According to the Agierre's manufacturing philosophy the VacuumConveyors with C.I.P. system can be customized according to the specific needs of each Customer.

Sistemi di trasporto in versione FILTERVAC Vacuum Transfer System in FILTERVAC version

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The Vacuum Conveyor type FILTERVAC is used to charge machines like Mixer, Reactor and Bicone where it's not possible to install the Vacuum Conveyor directly on the machine. The main functioning is based on the product suction by vacuum created inside the machine.

The Pump installed on the FILTERVAC generates the vacuum in the Machine that allows the continuous product suction from a drum, bin or other machines.

During the charging phase the FILTERVAC makes also the filter cleaning. At the end of the suction process is made the so called "purge phase" consisting in the powder recovery that can be done manually or automatically through the suction hopper installed under the discharge point of the FILTERVAC.

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Through the Ferrule of Vapour Collector the Vacuum Conveyor can be installed directly on reactors and/or dissolvers even with both solvents and vapour at 100°, thanks to a structure with double cavity wall which creates a barrier of laminar flow. The double cavity wall traps the vapour of reactor which goes up to the inlet connection of reactor. This system avoids any product accumulation on the discharge valve of the Vacuum Conveyor installed on the reactor, simplifying the discharging phase of product.

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The Lump Breaker is manufactured in accordance with the Machinery Directive 2006/42/EC and upon request it can be manufactured for installation in ATEX zones.

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The Lump Breaker model AGR-BREAKER allows to crush the lumps that could form in the materials during the production process or the handling of powders. The Lump Breaker Systems are particularly suitable for hygroscopic materials or those that tend to clog, and thanks to their versatility can be used in all industrial sectors.

It is characterized by a self-cleaning design with very small dimensions and can also be mounted within existing systems or systems for emptying (bags, drums and Big-Bag), or usually located immediately above vessel, filter dryers, centrifuges, mixers and generally on the process machine.

The breaking of lumps takes place through the rotation of special blades within a grid of contrast moved by means of a gearmotor.

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The lump breaker aspiration lance is designed to work in vacuum and to be connected with our Vacuum Conveyor. It finds the main application in the chemical and pharmaceutical industries for the handling of humid powders with big lumps that the standard suction lance can't transport.

The lump breaker lance can be manufactured also with a fluidisation air/nitrogen system which guarantee the right air/nitrogen proportion for the fluid dynamics balance.

According to the Agierre's manufacturing philosophy the lump breaker lance is manufactured in compliance with the "cGMP guidelines", UNI EN ISO and safety norms.

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Vacuum Transfer System PED (Pressure Rated)

AGIERRE[®]
Engineering & Innovative Solutions

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Vacuum Transfer System Pressure Rated

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All Vacuum Conveyors PED are manufactured according to the "cGMP guidelines"; the UNI EN ISO norms, ATEX Directive, Machinery Directive and PED Directive 97/23/EC.

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The PED Conveyors are manufactured for the need to install the Vacuum Conveyor System on reactors PED classified (according to the Directive 97/23/EC).

The main specifications are:

- Material: AISI 316L (Stainless Steel);
- Thickness of plate: 30/10;
- Internal and external finish: mirror polished;
- Connection between body and lid: flange, wing nut or clamp;
- Functioning: pneumatic, electropneumatic or with PLC
- Pressure: -1+3 bargs / -1+6 bargs

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
Also the PED Conveyors can be manufactured with C.I.P. (Cleaning in place) and S.I.P. (Steam in place) Systems. This type of Vacuum Conveyor is named GMA.

- | | |
|---|--|
| <p>1 Connessione Pompa e ingresso prodotto : Valvola a farfalla
<i>Connection Vacuum Pump and inlet connection : Butterfly valve</i></p> <p>2 Connessione aria/azoto in pressione: Valvola a Membrana/Farfalla
<i>Connection air/nitrogen: Diaphragm/Butterfly Valve</i></p> <p>3 Valvola di Sicurezza
<i>Safety Valve</i></p> | <p>4 Filtro Acciaio Inox
<i>Stainless Steel Filter</i></p> <p>5 Filtro in Acciaio inox con sistema di lavaggio per C.I.P. e S.I.P.
<i>Stainless Steel Filter with the system of cleaning for C.I.P. and S.I.P.</i></p> <p>6 Valvola a Farfalla sanitaria per ambienti sterili
<i>Sanitary Butterfly Valve for Sterile Areas</i></p> |
|---|--|

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According to the AGIERRE's projecting and manufacturing philosophy, also this type of Vacuum Conveyor (PED) can be customized according to the needs of the Customer.

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Vacuum Transfer System in HOPPERVAC version

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The Vacuum Conveyor type HOPPERVAC is mainly used when it's necessary to transport the product into a tank and/or container using the HOPPERVAC as a storage unit to supply another machine. In this case the HOPPERVAC is used not only as a powder transfer system but also as a container of product accumulation.

The functioning is the same of standard Vacuum Conveyor based on the product suction by the vacuum pump.

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
The HOPPERVAC like all other AGIERRE's Vacuum Transfer Systems can be manufactured with the C.I.P. system (Cleaning in Place) which allows the complete cleaning of the filter and the internal walls of Vacuum Conveyor.

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HOPPERVAC ro
charge mill and
dosing system

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According to the AGIERRE's projecting philosophy, the HOPPERVAC can be customized according to the needs of customer, always in compliance with the "cGMP guidelines", the UNI EN ISO norm, the Machine and ATEX Directives and with the safety norm.

Applicazioni Conveyor in Vacuum *Vacuum Conveyor Installations*